

Glass Nebulizers

Conikal

2mL/min uptake

[ARG-07-UC2](#)**SeaSpray**

2mL/min uptake

[ARG-07-USS2](#)**Slurry**

4mL/min uptake

[ARG-07-US6](#)**MicroMist**

600µL/min uptake

[ARG-07-UM06](#)

200µL/min uptake

[ARG-07-UM02](#)

100µL/min uptake

[ARG-07-UM01](#)

50µL/min uptake

[ARG-07-UM005](#)

Nebulizers suitable for use with HF

OpalMist

2mL/min uptake

[ARG-07-PFA2](#)

600µL/min uptake

[ARG-07-PFA06](#)

200µL/min uptake

[ARG-07-PFA02](#)

100µL/min uptake

[ARG-07-PFA01](#)

50µL/min uptake

[ARG-07-PFA005](#)**DuraMist**

1mL/min uptake

[ARG-07-DM1](#)

400µL/min uptake

[ARG-07-DM04](#)**VeeSpray - Ceramic**

Standard

[ARG-07-CV6](#)

IsoMist

IsoMist Programmable Temperature Spray Chamber

With glass spray chamber

[KT-1027](#)

With quartz spray chamber

[KT-1027Q](#)


With PFA spray chamber

[KT-1027P](#)

Spray Chambers

Twister

50mL Cyclonic

[20-809-2951HE](#)**Straight Quartz Adaptor**[21-809-3773](#)**Straight Quartz Adaptor with Aux Port**[21-809-3804](#)**Tracey TFE**

50mL Cyclonic

[20-809-3717](#)

Suitable for use with HF.

Requires Transfer Tube [21-809-3726](#).**Straight HF Resistant Adaptor**[21-809-3727](#)

Torches

Quartz Torch[31-808-3405](#)**Injectors**

Quartz Injector 1.2mm

[31-808-3410](#)

Quartz Injector 1.8mm

[31-808-3391](#)

Quartz Injector 2.0mm

[31-808-3411](#)

Quartz Injector 3.0mm

[31-808-3412](#)

Alumina Injector 1.8mm

[31-808-3409](#)

Alumina Injector 2.5mm

[31-808-3390](#)